

# 1 A

## **Cannula Cricothyroidotomy Kit**

## **Cannula Cricothyroidotomy Kit**

### **Contents**

14G cannula  
5 ml Luer lock syringe  
10 ml saline

### **Instructions**

Fill syringe with 5ml saline  
Remove cap from blunt end of cannula  
Connect saline-filled syringe to cannula

# 1 B

## Oxygenation Device

## Oxygenation Device

### Contents

O2 Tubing & Connector  
Swivel Adaptor  
Extension Tubing  
Luer Lock

### Instructions

Connect extension tubing to Cannula  
Insufflate Oxygen at HIGH FLOW  
Use finger to occlude swivel connector Y  
Oxygenate ON : OFF 3 secs : SpO2 drop ~5%

# 1 C

## Mini Trach Kit

## Mini Trach Kit

### Instructions

### WAIT BEFORE OPENING

When ready...,

- Have guidewire ready
- Have scalpel ready
- Have Mini Trach kit ready

# 2

## **Scalpel - Bougie Kit**

## **Scalpel - Bougie Kit**

### **Contents**

Size 10 Scalpel  
Frova Bougie  
Rapi-fit Connectors  
Size 6.0 Endotracheal Tube

### **Instructions**

Open all kit and pass to Doctor in order above

# 3

## Scalpel - Finger Kit

## Scalpel - Finger Kit

### Contents

Size 10 Scalpel  
Green Gauze  
14G Cannula  
5 ml Luer lock syringe  
10 ml Saline

### Instructions

Have scalpel ready  
Rest as per 1a kit

