

Before RSI performed



Equipment



Team brief



After RSI performed

- ✓ Airway assistant
- ✓ M.I.L.S / RAMP
- ✓ Drug giver & vocalise time
- ✓ ApO<sub>2</sub> on high flow
- ✓ Consider cricoid pressure

- ✓ Fluids runs easily
- ✓ RSI drug doses
- ✓ BVM with ETCO<sub>2</sub>
- ✓ OPA & NPA
- ✓ Laryngoscopes x 2
- ✓ Blade size
- ✓ Suction
- ✓ Bougie or stylet
- ✓ ETT & cuff tested
- ✓ Syringe for cuff
- ✓ Tube tie
- ✓ LMA (inc iLMA)
- ✓ Surgical airway kit

- ✓ Optimise patient position
- ✓ Maximise 360° access
- ✓ Pre-oxygenation
- ✓ Nasal prongs
- ✓ O<sub>2</sub> cylinder > half full & backup
- ✓ Baseline obs
- ✓ BP set at 2 minute cycles
- ✓ ETCO<sub>2</sub> connected to monitor, working

- ✓ Waveform ETCO<sub>2</sub> confirmed
- ✓ Check obs
- ✓ Tube secured & position noted
- ✓ Ventilation
- ✓ Sedation
- ✓ Paralysis
- ✓ Analgesia
- ✓ Oro/naso gastric tube
- ✓ Temperature probe

**CHALLENGE & RESPONSE**